### **METHOD STATEMENT**



#### **Determinand:**

Bromide, Sulphate, Chloride, Nitrite and Nitrate.

#### **Matrix:**

Surface waters, groundwaters, prepared soil leachates, trade effluents, treated sewage effluents, untreated sewage and landfill leachates for Sulphate, Chloride and Bromide. Then for Nitrite, Nitrate and Sulphate process waters and recreational waters.

### **Principle of Method:**

The samples are analysed using an ion chromatography system, fitted with a conductimetric detector. Automatic calculation is performed by the system's data analysis software.

#### **Sampling and Sample Preparation**

Samples are received in 1 litre PET bottles for sulphate, chloride, nitrate, and nitrite analysis. Subsamples are taken from the 1L PET bottle into a 10ml plastic sample tube. 100ml amber bottles (STL 061) are received for bromide analysis. All samples to be stored at 3±2°C

Samples are stable for times stated below, from sampling: Bromide 28 Days (ISO 5667:3)
Sulphate 30 Days (In-House Data)
Chloride 28 Days (ISO 5667:3)
Nitrite/TON 18 Days (In House Data)
Nitrate 30 days (In House Data)

#### Interferences:

Any substance with the same retention time as the analytes on the ion chromatograph may cause interference.

#### **Performance of Method:**

Range of Application without dilution: - Sulphate 0.20 - 400 mg/l

Bromide 0.20 - 50 mg/l Chloride 0.20 - 400 mg/l Nitrite 1.00 - 200mg/l 0.20 - 400 mg/l Nitrate Nitrite as N 0.30 - 60.87 mg/l 0.05 - 90.32 mg/l Nitrate as N 1.50-300 mg/l Nitrite as NaNO<sub>2</sub> 0.35 - 151.19 mg/l TON as NO<sub>3</sub>

# **METHOD STATEMENT**



### ICS2000:

Determinand	LOD	MRL	MRL Low Standard		High Standard		
Determinand	(mg/l)	(mg/l)	% RSD	% Bias	% RSD	% Bias	
Bromide	0.0188	0.20	2.70	1.36	1.69	0.90	
Sulphate	0.142	0.20	1.64	3.64	0.76	4.47	
Chloride	0.059	0.20	1.04	1.49	0.91	2.84	
Nitrite	0.0712	1.00	1.47	-0.31	1.72	-1.17	
Nitrate	0.0298	0.20	2.17	0.23	1.25	2.48	

Matrix	Sulphate		Bromide		Chloride	
IVIdUIX	% Recovery	% RSD	% Recovery	% RSD	% Recovery	% RSD
Groundwater	98.44	0.79	97.08	1.97	99.33	0.80
Surface Water	97.15	0.80	99.94	1.15	98.89	0.90
Final Effluent	98.03	0.71	99.29	1.63	98.78	0.99
Soil Leachate	98.87	0.79	92.70	2.11	99.30	0.79
Landfill Leachate	97.84	0.77	99.09	1.32	99.26	0.91
Untreated Sewage	96.63	0.69	99.56	1.63	97.21	0.82
Trade Effluent	100.14	0.94	99.44	1.30	98.78	20.99

Matrix	Nitrite		Nitrate		Sulphate	
ividuix	% Recovery	% RSD	% Recovery	% RSD	% Recovery	% RSD
Clean Process	101.60	2.22	103.29	1.12	102.53	1.26
Dirty Process	103.66	1.33	100.78	0.95	101.03	0.87
Recreational	99.69	1.16	101.43	0.86	100.84	0.71

Integrion:

LOD	MRL	Low Sta	andard	High Standard			
Determinand	(mg/l)	(mg/l)	(mg/l)	% RSD	% Bias	% RSD	% Bias
Bromide	0.0344	0.20	1.62	10.14	2.38	3.61	
Sulphate	0.1421	0.20	1.48	8.60	2.37	0.67	
Chloride	0.0706	0.20	1.50	3.27	2.40	0.13	
Nitrite	0.0466	1.00	1.34	2.84	2.80	1.93	
Nitrate	0.1236	0.20	1.55	2.18	2.37	1.50	

## **METHOD STATEMENT**



Matrix	Sulphate		Bromide		Chloride	
IVIALITX	%Recovery	%RSD	%Recovery	%RSD	%Recovery	%RSD
Surface Water	98.19	2.01	100.89	1.93	92.60	1.94
Trade Effluent	97.77	1.77	102.16	1.67	99.64	1.71
Final Effluent	98.19	1.97	101.42	1.88	90.87	1.94
Untreated Sewage	98.50	1.68	97.81	1.57	91.20	1.64
Ground water	99.58	1.88	100.60	2.06	94.83	1.97
Prepared Leachate	99.99	2.04	102.26	1.98	100.18	1.88
Landfill Leachate	98.01	1.93	96.82	2.21	95.28	2.05

Matrix	Nitrite		Nitrate		Sulphate	
IVIALIX	%Recovery	%RSD	%Recovery	%RSD	%Recovery	%RSD
Clean Process	100.57	2.57	100.26	2.25	100.98	2.36
Dirty Process	95.38	3.00	98.85	2.20	96.03	1.98
Recreational	100.11	2.70	100.26	2.67	102.91	2.37

### **Uncertainty of Measurement:**

The reported uncertainty is an expanded uncertainty calculated using a coverage factor of 2, which gives a level of confidence of approximately 95%.

Determinand	Uncertainty of Measurement (%)
Bromide	9.26
Sulphate	8.95
Chloride	13.68
Nitrite	9.08
Nitrate	6.54

### **References:**

The determination of Anions and Cations, Transition Metals, Other Complex Ions and Organic Acids and bases in Water by chromatography 1990 ISBN 0-11-752331-3.

ISO 5667-3:2024- Water quality Sampling Part 3: Preservation and handling of water samples.