## **METHOD STATEMENT**



#### **Determinand:**

Dissolved methane

#### **Matrix:**

Groundwaters, water extracts and leachates.

## **Principle of Method:**

A known amount of water is sealed in a headspace vial. The sample is then heated and agitated. This drives the dissolved methane into the 'headspace' above the aqueous layer. A portion of the headspace can then be sampled and analysed by GC using Flame Ionisation Detection.

# **Sampling and Sample Preparation**

This test is carried out on the water sample as received from the client.

The sample should be stored in a 40ml glass vial in a cold store at  $3 \pm 2^{\circ}$ C until required.

#### Interferences:

Any compound, eluting at the same retention time as methane and giving a response on the flame ionisation detector has the potential to interfere with the result.

### **Performance of Method:**

Range of Application: 0.01 - 25 mg/l

Normal Reporting Limit: 0.1mg/l

Determinand	LOD	Low St	andard	High Standard		
Determinand	(mg/l)	RSD %	Bias %	RSD %	Bias %	
Methane	0.00211	5.67	2.35	5.34	-1.10	

	Leachate			Surface Water			Ground Water					
Determinand	Low Spike H		High	n Spike Low		Spike	High Spike		Low Spike		High Spike	
	RSD %	Rec %	RSD %	Rec %	RSD %	Rec %	RSD %	Rec %	RSD %	Rec %	RSD %	Rec %
Methane	3.33	99.9	3.85	97.5	6.10	95.9	5.41	96.3	4.33	96.3	5.68	95.9

## **Uncertainty of Measurement:**

The reported uncertainty is an expanded uncertainty calculated using a coverage factor of 2, which gives a level of confidence of approximately 95%.

Determinand	Uncertainty of Measurement %
Methane	14.22

## **References:**

The determination of Methane and other Hydrocarbon Gases in Water 1988. ISBN 0117521280.

Method 62 Last Updated: 04.03.22 Page 1 of 1