METHOD STATEMENT



Determinand:

HBCDD

Matrix:

Surface, Ground and Saline waters

Principle of Method:

Approximately 400mls of sample is treated with sodium tetraethylborate (STEB) to ethylate TBT cation and is then extracted with 20mls of iso-hexane. The extract is analysed with a portion of the hexane extract exchanged into 70% methanol by LC-MSMS.

Sampling and Sample Preparation:

Samples should be taken in a 500ml amber or green glass bottle with 25mls of methanol as a preservative.

Samples are stored at $3 \pm 2^{\circ}$ C and analysed with 14 days of sampling. (In-house stability data)

Interferences:

MS/MS is an extremely selective technique and interferences should only be encountered very rarely, however in theory, any compound which is extracted by the procedure, which has a chromatographic retention time similar to the compound of interest and which produces both precursor and product ions similar to that of the compounds in question, may interfere

Performance of Method:

Compound	LOD (ng/L)	Low St	andard	High Standard	
		Bias	RSD	Bias	RSD
alpha-HBCDD	0.0416	-4.30%	11.43%	4.70%	8.09%
beta-HBCDD	0.0465	-4.47%	6.25%	2.69%	3.08%
gamma-HBCDD	0.0279	-4.73%	4.44%	3.20%	1.94%
HBCDD - Total	0.1378	-3.86%	6.30%	4.14%	3.22%

Precision, Bias and Limit of Detection

Matrix Spike Recoveries

Compound	Surface Water		Ground Water		Saline Water	
	Recovery	RSD	Recovery	RSD	Recovery	RSD
alpha-HBCDD	109.85%	9.15%	109.30%	10.01%	112.48%	13.58%
beta-HBCDD	103.27%	6.08%	101.40%	4.61%	102.01%	7.18%
gamma-HBCDD	103.33%	2.31%	102.50%	3.25%	102.95%	2.54%
HBCDD - Total	105.67%	3.69%	104.59%	4.11%	106.13%	3.09%

METHOD STATEMENT



Uncertainty of Measurement

The reported uncertainty is an expanded uncertainty calculated using a coverage factor of 2, which gives a level of confidence of approximately 95%.

Compound	Relative UoM (%)	Minimum UoM (ug/L)
alpha-HBCDD	23.9	0.000027
beta-HBCDD	27.6	0.000027
gamma-HBCDD	14.4	0.000020
HBCDD - Total	13.4	0.00008

References:

EU Priority Substances Directive 2013 Directive 2013/39/EC