



## Method Summary

### Determination of BTEX and GRO in Waters

#### Scope and Range

This method describes a procedure for determining the concentration of Gasoline Range Organic Hydrocarbons (GRO), in the carbon number range C<sub>5</sub> to C<sub>12</sub>, including BTEX components namely Benzene, Toluene, Ethylbenzene and Xylene isomers (o-, m- and p- isomers), and MTBE (tert butyl methyl ether) & TAME (tert amyl methyl ether). Method applies to waters and leachates.

Analyte	LOD (ug/l)	Upper Limit	Treated Sewage Accreditation	Ground Water Accreditation	Surface Water Accreditation
MTBE	<3	4000	ISO 17025	ISO 17025	ISO 17025
TAME	<4	4000	N/A	N/A	ISO 17025
Benzene	<7	4000	ISO 17025	ISO 17025	ISO 17025
Toluene	<4	4000	ISO 17025	ISO 17025	ISO 17025
Ethylbenzene	<5	4000	ISO 17025	ISO 17025	ISO 17025
m,p-Xylene	<8	4000	ISO 17025	ISO 17025	ISO 17025
o-Xylene	<3	4000	ISO 17025	ISO 17025	ISO 17025
GRO >C5-C12	<50	72,000	ISO 17025	ISO 17025	ISO 17025
GRO >C6-C8	<10	72,000	ISO 17025	ISO 17025	N/A
GRO >C8-C10	<10	72,000	ISO 17025	ISO 17025	N/A
GRO >C10-C12	<10	72,000	ISO 17025	ISO 17025	N/A
<b>A range of unaccredited aliphatic and aromatic bands are also reportable.</b>					

#### References

none

#### Principle

Preparation and Extraction:

Sodium Chloride is added to a clean 10ml headspace vial. Sample is added to the vial and surrogate standard is spiked beneath the surface of the water. The cap is then crimped in place.

Analysis:

Analysis is undertaken by gas chromatography with flame ionisation detection (GC-FID), e.g. - Agilent 6850 gas chromatograph with Gerstel or CTC headspace auto sampler and FID detector.

#### Interferences

Flame ionisation is a non-specific means of detection, therefore any substances either petroleum or non-petroleum in origin, that co-elute with any of the target compounds, will interfere with the determination e.g. solvents.