

METHOD STATEMENT

Determinand:

Determination of acid extractable fluoride

Matrix:

Sample Type: soil and sludges

Principle of Method:

A sample of soil (dried and ground) or sewage sludge (liquid) is extracted using 105°C dilute sulphuric acid. A citrate buffer is used to adjust the final pH to 5.4 ± 0.2 , release any complexed fluoride ions and to maintain the provision of constant high ionic strength. A fluoride selective electrode is used to measure the activity of fluoride in the prepared extract by means of standard addition technique, since the samples analysed contain substances likely to interfere.

Interferences:

Fluoride ions readily form complexes with polyvalent cations such as Al^{3+} , Ca^{2+} , Fe^{3+} , Mg^{2+} and PO_4^{3-} making them undetectable by ion selective electrode. Citrate buffer contains complexing agents, which release fluoride ions from metal-fluoride complexes that form. Under acidic conditions the formation of HF reduces the concentration of fluoride ions. The citrate buffer is designed to avoid these interferences by maintaining the pH of the sample between 5.2 and 5.6

Performance of Method:

Range of Application:

Reporting limit: 10mg/kg as F

Limit of Detection:

2.7mg/kg as F

Bias and Recoveries of Compounds:

Sample type	Mean sample result (mg/kg)	Bias (%)	RSD (%)	Spike recovery (%)
Soil	60.45	-	4.98	88.55
Liquid sludge	106.30	-	5.43	96.05

Sample type	RSD (%)	Prediction interval, mg/kg	Mean result, mg/kg
Sandy loam CRM	7.71	45.5 - 193	148.20
Clay CRM	7.51	77.9 - 412	291.80

References:

Fluoride in waters, effluents, sludges, plants and soils 1982; Methods for the Examination of Waters and Associated Materials; HMSO 1983 (ISBN 011 751662 7).

