

METHOD STATEMENT



Determinand:

Escherichia coli

Matrix:

Sample Type: sludges and soils

Principle of Method:

The sample is homogenised and serially diluted with maximum recovery diluent (MRD) and filtered through a membrane. The membrane filter is placed on an MLGA plate and *E.coli* are enumerated after incubation at $30 \pm 1.0^{\circ}\text{C}$ for 4 ± 0.5 hours, followed by $44 \pm 0.5^{\circ}\text{C}$ for 14 hours.

The MLGA contains lactose, an acidity indicator-phenol red and the chromagenic substrate BCIG (5-bromo-4-chloro-3-indolyl- β -glucuronide) either as a sodium or cyclohexylammonium salt. When hydrolysed BCIG indicates the presence of β -glucuronidase. Colonies that are β -glucuronidase positive are regarded as *E.coli*.

Sampling and Sample Preparation:

Samples should be taken in plastic containers. A minimum of 100g of sample should be provided to the laboratory.

Once taken, microbiological samples should be transferred immediately to dark storage conditions and kept at a temperature between $2 - 8^{\circ}\text{C}$ for transport to the laboratory. If samples are not analysed immediately on receipt in the laboratory, they should be kept at a temperature between $2 - 8^{\circ}\text{C}$ until analysis commences.

See the Microbiology of Sewage Sludge 2003, Part 2 - Practices and procedures for sampling and sample preparation (15.2) for more specific instruction on sampling.

Interferences

Interferences may include high levels of competing organisms or sample turbidity / debris.

Limitations - this method excludes a proportion of strains of *Escherichia coli* that are unable to grow at 44°C , do not express β -glucuronidase activity on primary isolation or are β -glucuronidase negative.

This method is not suitable for sludge samples that have been lime-treated or where enhanced microbial reduction is expected. These samples should be examined using a MPN technique.

The maximum number of colonies that should be counted on a filter membrane is 100. Anything reported higher than this is an estimation.

Performance of Method:

Range of Application:

The number of *Escherichia coli* present in sludge is expressed on a wet or dry weight basis. Typically results are reported as cfu of *Escherichia coli* per g of dried sludge

Limit of Detection:

The Limit of Detection has not yet been determined for this method

Uncertainty of measurement:

The Uncertainty of Measurement has not yet been determined for this method

References:

The Microbiology of Sewage Sludge (2003) - Part 1 - An overview of the treatment and use in agriculture of sewage sludge in relation to its impact on the environment and public health.

METHOD STATEMENT



Environment Agency - Methods for the Examination of Waters and Associated Materials. (Downloadable pdf only - no ISBN assigned).

The Microbiology of Sewage Sludge (2003) - Part 2 - Practices and procedures for sampling and sample preparation. Environment Agency - Methods for the Examination of Waters and Associated Materials. (Downloadable pdf only - no ISBN assigned).

The Microbiology of Sewage Sludge (2003) - Part 3 - Methods for the isolation and enumeration of *Escherichia coli*, including verocytotoxigenic *Escherichia coli*. Environment Agency - Methods for the Examination of Waters and Associated Materials. (Downloadable pdf only - no ISBN assigned).

The Microbiology of Sewage Sludge (2004) - Part 4 - Methods for the detection, isolation and enumeration of *Salmonellae*. Environment Agency - Methods for the Examination of Waters and Associated Materials. (Downloadable pdf only - no ISBN assigned).